

An easy way  
to maintain the sterility  
of the culture for a  
long time...

# LAMBDA MINIFOR Fermentor-Bioreactor

## Easy Sterility Concept

### How is it possible?





By eliminating the critical points causing contamination in common fermentors and bioreactors.

LAMBDA MINIFOR offers:

- ✓ Complete physical closure of the vessel and the side necks
- ✓ LAMBDA double-seal PEEK tubing connector
- ✓ Permanent stoppers with multiple seals
- ✓ High quality Peristaltic Pumps for long term continuous runs

### Unique construction for easy sterility



LAMBDA's unique sterility construction assures you with perfect and easy sterility maintenance for each and every run...



at LOWEST COST